

## 兔抗 FAM91A1 多克隆抗体

中文名称：兔抗 FAM91A1 多克隆抗体

英文名称：Anti-FAM91A1 rabbit polyclonal antibody

别名：family with sequence similarity 91 member A1

相关类别：一抗

储存：冷冻（-20℃）

抗原：FAM91A1

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	FAM91A1 (family with sequence similarity 91, member A1) is a 94 kDa phospho-protein that is not yet characterized.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	FAM91A1
<b>Immunogen:</b>	Synthetic peptide of human FAM91A1
<b>Full name:</b>	family with sequence similarity 91 member A1
<b>SwissProt:</b>	Q658Y4
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer and Human cervical cancer
<b>IHC Recommend dilution:</b>	30-150



