

兔抗 FAM98A 多克隆抗体

- 中文名称：兔抗 FAM98A 多克隆抗体
- 英文名称：Anti-FAM98A rabbit polyclonal antibody
- 别名：family with sequence similarity 98 member A
- 相关类别：一抗
- 抗原：FAM98A
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 反应种属：Human, Mouse, Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	FAM98A belongs to the FAM98 family. The function of FAM98A remains unknown.
Applications:	ELISA, WB, IHC
Name of antibody:	FAM98A
Immunogen:	Synthetic peptide of human FAM98A
Full name:	family with sequence similarity 98 member A
SwissProt:	Q8NCA5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	30-150
WB Predicted band size:	55 kDa

WB Positive control:	Human breast cancer tissue , Human liver cancer tissue and Human lung adenocarcinoma tissue lysates
WB Recommended dilution:	200-1000



