

兔抗 FAM110D 多克隆抗体

中文名称：兔抗 FAM110D 多克隆抗体

英文名称：Anti-FAM110D rabbit polyclonal antibody

别名：family with sequence similarity 110 member D; GRRP1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FAM110D

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, IHC
Name of antibody:	FAM110D
Immunogen:	Synthetic peptide of human FAM110D
Full name:	family with sequence similarity 110 member D
Synonyms:	GRRP1
SwissProt:	Q8TAY7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	30-150

