

## 兔抗 FAM129B 多克隆抗体

中文名称：兔抗 FAM129B 多克隆抗体

英文名称：Anti-FAM129B (Ab-712) rabbit polyclonal antibody

别名：OC58; MEG-3; C9orf88; MINERVA; bA356B19.6

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：FAM129B

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	a maternal imprinting gene. May possess growth suppressor activity. Two major alternatively spliced isoforms have been described. forms. One form, Meg3-proximal (Meg3p), contains exons 1-3. The second form, Meg3-distal (Meg3d) did not contain exons 1-3 and was present in oocytes and in 1- and 2-cell embryos
<b>Applications:</b>	WB, IF
<b>Name of antibody:</b>	FAM129B
<b>Immunogen:</b>	Synthetic peptide of human FAM129B (Ab-712)
<b>Full name:</b>	family with sequence similarity 129, member B
<b>Synonyms :</b>	OC58; MEG-3; C9orf88; MINERVA; bA356B19.6

<b>SwissProt:</b>	Q96TA1
<b>WB Predicted band size:</b>	84 kDa
<b>WB Positive control:</b>	3T3 and DU145 cells
<b>WB Recommended dilution:</b>	500-1000
<b>IF Positive control:</b>	Hela cells
<b>IF Recommended dilution</b>	100-200



