

兔抗 FAM177B 多克隆抗体

中文名称：兔抗 FAM177B 多克隆抗体

英文名称：Anti-FAM177B rabbit polyclonal antibody

别名：family with sequence similarity 177 member B

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FAM177B


反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, IHC
Name of antibody:	FAM177B
Immunogen:	Synthetic peptide of human FAM177B
Full name:	family with sequence similarity 177 member B
SwissProt:	A6PVY3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human tonsil
IHC Recommend dilution:	50-200



mlbio 梅联生物
Good elisakit producers

