

兔抗 FAM200A 多克隆抗体

- 中文名称：兔抗 FAM200A 多克隆抗体
- 英文名称：Anti-FAM200A rabbit polyclonal antibody
- 别名：family with sequence similarity 200 member A; C7orf38
- 抗原：FAM200A
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein of unknown function. The protein is weakly similar to transposase-like proteins in human and mouse.
Applications:	ELISA, IHC
Name of antibody:	FAM200A
Immunogen:	Synthetic peptide of human FAM200A
Full name:	family with sequence similarity 200 member A
Synonyms:	C7orf38
SwissProt:	Q8TCP9
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human thyroid cancer

IHC Recommend dilution:	50-300
-------------------------	--------



