

兔抗 FAM234A 多克隆抗体

中文名称: 兔抗 FAM234A 多克隆抗体

英文名称: Anti-FAM234A rabbit polyclonal antibody

别名: gs19; ITFG3; C16orf9

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: FAM234A

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	FAM234A, also known as ITFG3, which is a 552 amino acid single-pass transmembrane protein that contains FG-GAP repeats, a motif commonly found in Integrin proteins. There are two isoforms of ITFG3 that are produced as a result of alternative splicing events.
Applications:	ELISA, WB, IHC
Name of antibody:	FAM234A
Immunogen:	Synthetic peptide of human FAM234A
Full name:	family with sequence similarity 234 member A
Synonyms:	gs19; ITFG3; C16orf9
SwissProt:	Q9H0X4
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human brain
IHC Recommend dilution:	20-100
WB Predicted band size:	60 kDa
WB Positive control:	Rat brain tissue lysate
WB Recommended dilution:	200-1000



