

## 兔抗 FANCF 多克隆抗体

中文名称：兔抗 FANCF 多克隆抗体

英文名称： Anti-FANCF rabbit polyclonal antibody

别名： FAF

相关类别： 一抗

抗原： FANCF

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

The Fanconi anemia complementation group (FANC) currently includes FANCA, FANCB, FANCC, FANCD1 (also called BRCA2), FANCD2, FANCE, FANCF, FANCG, FANCI, FANCI (also called BRIP1), FANCL, FANCM and FANCN (also called PALB2). The previously defined group FANCH is the same as FANCA. Fanconi anemia is a genetically heterogeneous recessive disorder characterized by cytogenetic instability, hypersensitivity to DNA crosslinking agents, increased chromosomal breakage, and defective DNA repair. The members of the Fanconi anemia complementation group do not share sequence similarity; they are related by their as

	sembly into a common nuclear protein complex. This gene encodes the protein for complementation group F.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	FANCF
<b>Immunogen:</b>	Fusion protein of human FANCF
<b>Full name:</b>	Fanconi anemia, complementation group F
<b>Synonyms :</b>	FAF
<b>SwissProt:</b>	Q9NPI8
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human cervical cancer and Human brain
<b>IHC Recommend dilution:</b>	100-300

