

## XPO1 抗原（重组蛋白）

中文名称：XPO1 抗原（重组蛋白）

英文名称：XPO1 Antigen (Recombinant Protein)

别名：emb; CRM1; exp1; CRM-1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to C terminal 250 amino acids of human XPO1

技术规格：

<b>Full name:</b>	exportin 1
<b>Synonyms:</b>	emb; CRM1; exp1; CRM-1
<b>Swissprot:</b>	O14980
<b>Gene Accession:</b>	BC032847
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This cell-cycle-regulated gene encodes a protein that mediates leucine-rich nuclear export signal (NES)-dependent protein transport. The protein specifically inhibits the nuclear export of Rev and U snRNAs. It is involved in the control of several cellular processes by controlling the localization of cyclin B, MPAK, and MAPKAP kinase 2. This protein also regulates NFAT and AP-1.