

## 兔抗 FAP 多克隆抗体

中文名称：兔抗 FAP 多克隆抗体

英文名称： Anti-FAP rabbit polyclonal antibody

别名： FAPA; DPPIV

相关类别： 一抗

抗原： FAP

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a homodimeric integral membrane gelatinase belonging to the serine protease family. It is selectively expressed in reactive stromal fibroblasts of epithelial cancers, granulation tissue of healing wounds, and malignant cells of bone and soft tissue sarcomas. This protein is thought to be involved in the control of fibroblast growth or epithelial-mesenchymal interactions during development, tissue repair, and epithelial carcinogenesis.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	FAP

<b>Immunogen:</b>	Fusion protein of human FAP
<b>Full name:</b>	fibroblast activation protein, alpha
<b>Synonyms :</b>	FAPA; DPPIV
<b>SwissProt:</b>	Q12884
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human brain and Human thyroid cancer
<b>IHC Recommend dilution:</b>	50-200

