

兔抗 FAR2 多克隆抗体

中文名称： 兔抗 FAR2 多克隆抗体

英文名称： Anti-FAR2 rabbit polyclonal antibody

别名： MLSTD1; SDR10E2

相关类别： 一抗

储存： 冷冻 (-20℃)

抗原： FAR2

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene belongs to the short chain dehydrogenase/reductase superfamily. It encodes a reductase enzyme involved in the first step of wax biosynthesis wherein fatty acids are converted to fatty alcohols. The encoded peroxisomal protein utilizes saturated fatty acids of 16 or 18 carbons as preferred substrates. Alternatively spliced transcript variants have been observed for this gene. Related pseudogenes have been identified on chromosomes 2, 14 and 22.
Applications:	ELISA, IHC
Name of antibody:	FAR2

Immunogen:	Fusion protein of human FAR2
Full name:	fatty acyl CoA reductase 2
Synonyms :	MLSTD1; SDR10E2
SwissProt:	Q96K12
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human ovarian cancer and human tonsil
IHC Recommend dilution:	50-200

