

兔抗 FARS2 多克隆抗体

中文名称：兔抗 FARS2 多克隆抗体

英文名称：Anti-FARS2 rabbit polyclonal antibody

别名：phenylalanyl-tRNA synthetase 2, mitochondrial; FARS1; PheRS; SPG77; COXPD14; HSPC320; mtPheRS

相关类别：一抗

储存：冷冻 (-20℃)

抗原：FARS2

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein that transfers phenylalanine to its cognate tRNA. This protein localizes to the mitochondrion and plays a role in mitochondrial protein translation. Mutations in this gene can cause combined oxidative phosphorylation deficiency 14 (Alpers encephalopathy). Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB, IHC
Name of antibody:	FARS2
Immunogen:	Fusion protein of human FARS2

Full name:	phenylalanyl-tRNA synthetase 2, mitochondrial
Synonyms:	FARS1; PheRS; SPG77; COXPD14; HSPC320; mtPheRS
SwissProt:	O95363
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	150-300
WB Predicted band size:	52 kDa
WB Positive control:	293T and Hela cell lysates
WB Recommended dilution:	500-2000



