

兔抗 FASN 多克隆抗体

中文名称: 兔抗 FASN 多克隆抗体

英文名称: Anti-FASN rabbit polyclonal antibody

别名: FAS; OA-519; SDR27X1

抗原: FASN

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The enzyme encoded by this gene is a multifunctional protein. Its main function is to catalyze the synthesis of palmitate from acetyl-CoA and malonyl-CoA, in the presence of NADPH, into long-chain saturated fatty acids. In some cancer cell lines, this protein has been found to be used with estrogen receptor-alpha (ER-alpha), in which the N-terminus of FAS is fused in-frame with the C-terminus of ER-alpha.
Applications:	ELISA, IHC
Name of antibody:	FASN
Immunogen:	Synthetic peptide of human FASN

Full name:	fatty acid synthase
Synonyms :	FAS; OA-519; SDR27X1
SwissProt:	P49327
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	25-100

