

兔抗 FASTKD2 多克隆抗体

中文名称：兔抗 FASTKD2 多克隆抗体

英文名称：Anti-FASTKD2 rabbit polyclonal antibody

别名：KIAA0971

相关类别：一抗

抗原：FASTKD2

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein that is localized in the mitochondrial inner compartment and that may play a role in mitochondrial apoptosis. Nonsense mutations have been reported to result in cytochrome c oxidase deficiency.
Applications:	ELISA, IHC
Name of antibody:	FASTKD2
Immunogen:	Fusion protein of human FASTKD2
Full name:	FAST kinase domains 2
Synonyms :	KIAA0971
SwissProt:	Q9NYY8

ELISA Recommended dilution:	2000-10000
IHC positive control:	Human cervical cancer and human esophagus cancer
IHC Recommend dilution:	100-300

