

兔抗 FASTKD5 多克隆抗体

中文名称：兔抗 FASTKD5 多克隆抗体

英文名称：Anti-FASTKD5 rabbit polyclonal antibody

别名：dj1187M17.5

相关类别：一抗

抗原：FASTKD5

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

Representing about 2% of human DNA, chromosome 20 consists of approximately 63 million bases and 600 genes. Chromosome 20 contains a region with numerous genes expressed in the epididymis which are thought important for seminal production and some viewed as potential targets for male contraception. The PRNP gene encoding the prion protein associated with spongiform encephalopathies, like Creutzfeldt-Jakob disease, is found on chromosome 20. Amyotrophic lateral sclerosis, spinal muscular atrophy, ring chromosome 20 epilepsy syndrome and Alagille syndrome

	me are also associated with chromosome 20.
Applications:	ELISA, IHC
Name of antibody:	FASTKD5
Immunogen:	Fusion protein of human FASTKD5
Full name:	FAST kinase domains 5
Synonyms :	dJ1187M17.5
SwissProt:	Q7L8L6
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer and human brain
IHC Recommend dilution:	50-200

