

兔抗 FAT3 多克隆抗体

中文名称：兔抗 FAT3 多克隆抗体

英文名称：Anti-FAT3 rabbit polyclonal antibody

别名：CDHF15; CDHR10

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FAT3

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May play a role in the interactions between neurites derived from specific subsets of neurons during development.
Applications:	ELISA, IHC
Name of antibody:	FAT3
Immunogen:	Synthetic peptide of human FAT3
Full name:	FAT atypical cadherin 3
Synonyms:	hFat3; CDHF15; CDHR10
SwissProt:	Q8TDW7
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human esophagus cancer and Human gastric cancer

IHC Recommend dilution: 25-100

