

ABHD10 抗原（重组蛋白）

中文名称：ABHD10 抗原（重组蛋白）

英文名称：ABHD10 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to a region derived from 1-306 amino acids of human ABHD10

技术规格：

Full name:	abhydrolase domain containing 10
Swissprot:	Q9NUJ1
Gene Accession:	BC014516
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a mitochondrially-localized enzyme that acts in liver cells as a hydrolase. The encoded protein removes glucuronide from mycophenolic acid acyl-glucuronide. There is a pseudogene for this gene on chromosome 6. Alternative splicing results in multiple transcript variants.