

## 兔抗 FBLN2 多克隆抗体

中文名称：兔抗 FBLN2 多克隆抗体

英文名称： Anti-FBLN2 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： FBLN2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes an extracellular matrix protein, which belongs to the fibulin family. This protein binds various extracellular ligands and calcium. It may play a role during organ development, in particular, during the differentiation of heart, skeletal and neuronal structures. Alternatively spliced transcript variants encoding different isoforms have been identified.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	FBLN2
<b>Immunogen:</b>	Synthetic peptide of human FBLN2
<b>Full name:</b>	fibulin 2
<b>SwissProt:</b>	P98095
<b>ELISA Recommended dilution:</b>	5000-10000

<b>WB Predicted band size:</b>	126 kDa
<b>WB Positive control:</b>	Hela , A549 and 293T cell
<b>WB Recommended dilution:</b>	200-1000

