

兔抗 FBP1 多克隆抗体

中文名称：兔抗 FBP1 多克隆抗体

英文名称：Anti-FBP1 rabbit polyclonal antibody

别名：FBP

相关类别：一抗

储存：冷冻（-20℃）

抗原：FBP1

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Fructose-1,6-bisphosphatase 1, a gluconeogenesis regulatory enzyme, catalyzes the hydrolysis of fructose 1,6-bisphosphate to fructose 6-phosphate and inorganic phosphate. Fructose-1,6-diphosphatase deficiency is associated with hypoglycemia and metabolic acidosis.
Applications:	ELISA, WB, IHC
Name of antibody:	FBP1
Immunogen:	Fusion protein of human FBP1
Full name:	fructose-1,6-bisphosphatase 1
Synonyms :	FBP

SwissProt:	P09467
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human prostate cancer and human breast cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	37 kDa
WB Positive control:	Mouse stomach and human fetal liver tissue, MCF-7 cells and mouse liver tissue
WB Recommended dilution:	1000-5000



