

兔抗 FBP2 多克隆抗体

中文名称：兔抗 FBP2 多克隆抗体

英文名称：Anti-FBP2 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

抗 原：FBP2

宿 主：Rabbit

反应种属：Human, Mouse, Rat

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a gluconeogenesis regulatory enzyme which catalyzes the hydrolysis of fructose 1,6-bisphosphate to fructose 6-phosphate and inorganic phosphate.
Applications:	ELISA, IHC
Name of antibody:	FBP2
Immunogen:	Fusion protein of human FBP2
Full name:	fructose-1,6-bisphosphatase 2
SwissProt:	O00757
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human brain
IHC Recommend dilution:	50-200

