

ADRM1 抗原（重组蛋白）

中文名称：ADRM1 抗原（重组蛋白）

英文名称：ADRM1 Antigen (Recombinant Protein)

别名：adhesion regulating molecule 1; ARM1; ARM-1; GP110

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 208-407 amino acids of human ADRM1

技术规格：

Full name:	adhesion regulating molecule 1
Synonyms:	ARM1; ARM-1; GP110
Swissprot:	Q16186
Gene Accession:	BC017245
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the adhesion regulating molecule 1 protein family. The encoded protein is a component of the proteasome where it acts as a ubiquitin receptor and recruits the deubiquitinating enzyme, ubiquitin carboxyl-terminal hydrolase L5. Increased levels of the encoded protein are associated with increased cell adhesion, which is likely an indirect effect of this intracellular protein. Dysregulation of this gene has been implicated in carcinogenesis. Al

ternative splicing results in multiple transcript variants.