

AK9 抗原（重组蛋白）

中文名称：AK9 抗原（重组蛋白）

英文名称：AK9 Antigen (Recombinant Protein)

别名：AK 9; AKD1; AKD2; C6orf199; C6orf224; dJ70A9.1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 222-421 amino acids of human AK9

技术规格：

Full name:	adenylate kinase 9
Synonyms:	AK 9; AKD1; AKD2; C6orf199; C6orf224; dJ70A9.1
Swissprot:	Q5TCS8
Gene Accession:	BC022031
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene catalyzes the interconversion of nucleosides, possessing both nucleoside monophosphate and diphosphate kinase activities. The encoded protein uses these interconversions to maintain nucleoside homeostasis.