

AP2B1 抗原（重组蛋白）

中文名称：AP2B1 抗原（重组蛋白）

英文名称：AP2B1 Antigen (Recombinant Protein)

别名：ADTB2; AP105B; CLAPB1; AP2-BETA

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 738-937 amino acids of human AP2B1

技术规格：

Full name:	adaptor related protein complex 2 beta 1 subunit
Synonyms:	ADTB2; AP105B; CLAPB1; AP2-BETA
Swissprot:	P63010
Gene Accession:	BC006201
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is one of two large chain components of the assembly protein complex 2, which serves to link clathrin to receptors in coated vesicles. The encoded protein is found on the cytoplasmic face of coated vesicles in the plasma membrane. Two transcript variants encoding different isoforms have been found for this gene.