

ARPC5L 抗原（重组蛋白）

中文名称：ARPC5L 抗原（重组蛋白）

英文名称：ARPC5L Antigen (Recombinant Protein)

别名：ARC16-2

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-153 amino acids of human ARPC5L

技术规格：

Full name:	actin related protein 2/3 complex subunit 5 like
Synonyms:	ARC16-2
Swissprot:	Q9BPX5
Gene Accession:	BC000018
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May function as component of the Arp2/3 complex which is involved in regulation of actin polymerization and together with an activating nucleation-promoting factor (NPF) mediates the formation of branched actin networks.