

## ASTE1 抗原（重组蛋白）

中文名称：ASTE1 抗原（重组蛋白）

英文名称：ASTE1 Antigen (Recombinant Protein)

别名：HT001

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 480-679 amino acids of human ASTE1

技术规格：

|                           |  |
|---------------------------|--|
| <b>Full name:</b>         | asteroid homolog 1                                     |
| <b>Synonyms:</b>          | HT001  |
| <b>Swissprot:</b>         | Q2TB18   |
| <b>Gene Accession:</b>    | BC030274   |
| <b>Purity:</b>            | >85%, as determined by Coomassie blue stained SDS-PAGE |
| <b>Expression system:</b> | Escherichia coli                                       |
| <b>Tags:</b>              | His tag C-Terminus, GST tag N-Terminus                 |
| <b>Background:</b>        | Possible role in EGF receptor signaling.               |