

## ASTN2 抗原（重组蛋白）

中文名称：ASTN2 抗原（重组蛋白）

英文名称：ASTN2 Antigen (Recombinant Protein)

别名：bA67K19.1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1140-1339 amino acids of human ASTN2

技术规格：

<b>Full name:</b>	astrotactin 2
<b>Synonyms:</b>	bA67K19.1
<b>Swissprot:</b>	O75129
<b>Gene Accession:</b>	BC029272
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a protein that is expressed in the brain and may function in neuronal migration, based on functional studies of the related astrotactin 1 gene in human and mouse. A deletion at this locus has been associated with schizophrenia. Multiple transcript variants encoding different proteins have been found for this locus.

**mlbio** 梅联生物  
Good elisakit producers