

C19orf57 抗原（重组蛋白）

中文名称：C19orf57 抗原（重组蛋白）

英文名称：C19orf57 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 471-658 amino acids of human C19orf57

技术规格：

Full name:	chromosome 19 open reading frame 57
Swissprot:	Q0VDD7
Gene Accession:	BC012945
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The function of this protein remains unknown.