

CDC20 抗原（重组蛋白）

中文名称：CDC20 抗原（重组蛋白）

英文名称：CDC20 Antigen (Recombinant Protein)

别名：CDC20A; p55CDC; bA276H19.3

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 300-499 amino acids of human CDC20

技术规格：

Full name:	cell division cycle 20
Synonyms:	CDC20A; p55CDC; bA276H19.3
Swissprot:	Q12834
Gene Accession:	BC000624
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	CDC20 appears to act as a regulatory protein interacting with several other proteins at multiple points in the cell cycle. It is required for two microtubule-dependent processes, nuclear movement prior to anaphase and chromosome separation.