

## CFAP47 抗原（重组蛋白）

中文名称：CFAP47 抗原（重组蛋白）

英文名称：CFAP47 Antigen (Recombinant Protein)

别名：cilia and flagella associated protein 47; CHDC2; CXorf22; CXorf30; CXorf59

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-200 amino acids of human CFAP47

技术规格：

<b>Full name:</b>	cilia and flagella associated protein 47
<b>Synonyms:</b>	CHDC2; CXorf22; CXorf30; CXorf59
<b>Swissprot:</b>	Q6ZTR5
<b>Gene Accession:</b>	BC027936
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The function of this protein remains unknown.