

CNPY3 抗原（重组蛋白）

中文名称：CNPY3 抗原（重组蛋白）

英文名称：CNPY3 Antigen (Recombinant Protein)

别名：canopy FGF signaling regulator 3; CAG4A; ERDA5; TNRC5; PRAT4A

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 31-278 amino acids of human CNPY3

技术规格：

Full name:	canopy FGF signaling regulator 3
Synonyms:	CAG4A; ERDA5; TNRC5; PRAT4A
Swissprot:	Q9BT09
Gene Accession:	BC004423
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein that binds members of the toll-like receptor protein family and functions as a chaperone to aid in folding and export of these proteins. Alternative splicing results in multiple transcript variants. Naturally occurring readthrough transcription occurs between this locus and the downstream GNMT (glycine N-methyltransferase) gene and is represented with GeneID:107080644.