

兔抗 FDCSP 多克隆抗体

中文名称：兔抗 FDCSP 多克隆抗体

英文名称：Anti-FDCSP rabbit polyclonal antibody

别名：C4orf7; FDC-SP

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FDCSP

反应种属：Human


标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a small secreted protein that is expressed in follicular dendritic cells. This protein specifically binds to activated B cells, and functions as a regulator of antibody responses. It is also thought to contribute to tumor metastases by promoting cancer cell migration and invasion.
Applications:	ELISA, IHC
Name of antibody:	FDCSP
Immunogen:	Synthetic peptide of human FDCSP
Full name:	follicular dendritic cell secreted protein
Synonyms:	C4orf7; FDC-SP

SwissProt:	Q8NFU4
IHC positive control:	Human gastric cancer and human thyroid cancer
IHC Recommend dilution:	25-100



mlbio 梅联生物
Good elisakit producers