

## DCAF10 抗原（重组蛋白）

中文名称：DCAF10 抗原（重组蛋白）

英文名称：DCAF10 Antigen (Recombinant Protein)

别名：WDR32

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 389-559 amino acids of human DCAF10

技术规格：

<b>Full name:</b>	DDB1 and CUL4 associated factor 10
<b>Synonyms:</b>	WDR32
<b>Swissprot:</b>	Q5QP82
<b>Gene Accession:</b>	BC003520
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	May function as a substrate receptor for CUL4-DDB1 E3 ubiquitin-protein ligase complex.