

兔抗 FDFT1 多克隆抗体

中文名称：兔抗 FDFT1 多克隆抗体

英文名称：Anti-FDFT1 rabbit polyclonal antibody

别名：farnesyl-diphosphate farnesyltransferase 1; SS; SQS; DGPT; ERG9

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FDFT1

反应种属：Human, Mouse, Rat


标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a membrane-associated enzyme located at a branch point in the mevalonate pathway. The encoded protein is the first specific enzyme in cholesterol biosynthesis, catalyzing the dimerization of two molecules of farnesyl diphosphate in a two-step reaction to form squalene.
Applications:	ELISA, WB
Name of antibody:	FDFT1
Immunogen:	Fusion protein of human FDFT1
Full name:	farnesyl-diphosphate farnesyltransferase 1
Synonyms:	SS; SQS; DGPT; ERG9

SwissProt:	P37268
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	48 kDa
WB Positive control:	Human fetal brain tissue, HepG2 and HT29 cell lysates
WB Recommended dilution:	1000-5000



mlbio 梅联生物
Good elisakit producers