

## DENND6A 抗原（重组蛋白）

中文名称：DENND6A 抗原（重组蛋白）

英文名称：DENND6A Antigen (Recombinant Protein)

别名：AFI1A; FAM116A

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 409-608 amino acids of human DENND6A

技术规格：

<b>Full name:</b>	DENN domain containing 6A
<b>Synonyms:</b>	AFI1A; FAM116A
<b>Swissprot:</b>	Q8IWF6
<b>Gene Accession:</b>	BC040291
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Guanine nucleotide exchange factor (GEF) for RAB14. Component of a n endocytic recycling pathway that is required for the control of ADA M10 transport, shedding of N-cadherin/CDH2 by ADAM9 or ADAM10 and regulation of cell-cell junctions. Required for RAB14 recruitment to recycling endosomes.