

DNPEP 抗原（重组蛋白）

中文名称：DNPEP 抗原（重组蛋白）

英文名称：DNPEP Antigen (Recombinant Protein)

别名：aspartyl aminopeptidase; DAP; ASPEP

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 276-475 amino acids of human DNPEP

技术规格：

Full name:	aspartyl aminopeptidase
Synonyms:	DAP; ASPEP
Swissprot:	Q9ULA0
Gene Accession:	BC000653
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is an aminopeptidase which prefers acidic amino acids, and specifically favors aspartic acid over glutamic acid. It is thought to be a cytosolic protein involved in general metabolism of intracellular proteins. Several transcript variants encoding different isoforms have been found for this gene.