

## DPY19L3 抗原（重组蛋白）

中文名称：DPY19L3 抗原（重组蛋白）

英文名称：DPY19L3 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to a region derived from 1-102 amino acids of human DPY19L3

技术规格：

<b>Full name:</b>	dpy-19 like C-mannosyltransferase 3
<b>Swissprot:</b>	Q6ZPD9
<b>Gene Accession:</b>	BC029162
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Probable C-mannosyltransferase that mediates C-mannosylation of tryptophan residues on target proteins.