

## ELAC1 抗原（重组蛋白）

中文名称：ELAC1 抗原（重组蛋白）

英文名称：ELAC1 Antigen (Recombinant Protein)

别名：D29

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 164-363 amino acids of human ELAC1

技术规格：

<b>Full name:</b>	elaC ribonuclease Z 1
<b>Synonyms:</b>	D29
<b>Swissprot:</b>	Q9H777
<b>Gene Accession:</b>	BC014624
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Zinc phosphodiesterase, which displays some tRNA 3'-processing endonuclease activity. Probably involved in tRNA maturation, by removing a 3'-trailer from precursor tRNA.