

EMC1 抗原（重组蛋白）

中文名称：EMC1 抗原（重组蛋白）

英文名称：EMC1 Antigen (Recombinant Protein)

别名：CAVIPMR; KIAA0090

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 22-221 amino acids of human EMC1

技术规格：

Full name:	ER membrane protein complex subunit 1
Synonyms:	CAVIPMR; KIAA0090
Swissprot:	Q8N766
Gene Accession:	BC034589
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a single-pass type I transmembrane protein, which is a subunit of the endoplasmic reticulum membrane protein complex (EMC). Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene.