

兔抗 FGF21 多克隆抗体

中文名称：兔抗 FGF21 多克隆抗体

英文名称：Anti-FGF21 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

抗 原：FGF21

宿 主：Rabbit

反应种属：Human, Mouse

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities and are involved in a variety of biological processes including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. The function of this growth factor has not yet been determined.
Applications:	ELISA, IHC
Name of antibody:	FGF21
Immunogen:	Fusion protein of human FGF21
Full name:	fibroblast growth factor 21
SwissProt:	Q9NSA1

ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer and human brain
IHC Recommend dilution:	25-100

