

兔抗 FGF22 多克隆抗体

中文名称：兔抗 FGF22 多克隆抗体

英文名称：Anti-FGF22 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻 (-20℃)

宿 主：Rabbit

抗 原：FGF22

反应种属：Human, Mouse

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities and are involved in a variety of biological processes including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. The mouse homolog of this gene was found to be preferentially expressed in the inner root sheath of the hair follicle, which suggested a role in hair development.
Applications:	ELISA, IHC
Name of antibody:	FGF22
Immunogen:	Synthetic peptide of human FGF22
Full name:	fibroblast growth factor 22

SwissProt:	Q9HCT0
IHC positive control:	Human gastric cancer and human colon cancer
IHC Recommend dilution:	25-100

