

FDXR 抗原（重组蛋白）

中文名称：FDXR 抗原（重组蛋白）

英文名称：FDXR Antigen (Recombinant Protein)

别名：ADXR; ANOA

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 292-491 amino acids of human FDXR

技术规格：

Full name:	ferredoxin reductase
Synonyms:	ADXR; ANOA
Swissprot:	P22570
Gene Accession:	BC002960
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a mitochondrial flavoprotein that initiates electron transport for cytochromes P450 receiving electrons from NADPH. Multiple alternatively spliced transcript variants have been found for this gene.