

FUBP3 抗原（重组蛋白）

中文名称：FUBP3 抗原（重组蛋白）

英文名称：FUBP3 Antigen (Recombinant Protein)

别名：far upstream element binding protein 3; FBP3

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-261 amino acids of human FUBP3

技术规格：

Full name:	far upstream element binding protein 3
Synonyms:	FBP3
Swissprot:	Q96I24
Gene Accession:	BC007874
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression.