

兔抗 FHL3 多克隆抗体

中文名称：兔抗 FHL3 多克隆抗体

英文名称：Anti-FHL3 rabbit polyclonal antibody

别名：SLIM2

相关类别：一抗

抗原：FHL3

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is a member of a family of proteins containing a four-and-a-half LIM domain, which is a highly conserved double zinc finger motif. The encoded protein has been shown to interact with the cancer developmental regulators SMAD2, SMAD3, and SMAD4, the skeletal muscle myogenesis protein MyoD, and the high-affinity IgE beta chain regulator MZF-1. This protein may be involved in tumor suppression, repression of MyoD expression, and repression of IgE receptor expression. Two transcript variants encoding different isoforms have been found for this gene.

Applications:	ELISA, WB, IHC
Name of antibody:	FHL3
Immunogen:	Fusion protein of human FHL3
Full name:	four and a half LIM domains 3
Synonyms :	SLIM2
SwissProt:	Q13643
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human cervical cancer and human brain
IHC Recommend dilution:	50-200
WB Predicted band size:	31 kDa
WB Positive control:	Mouse skeletal muscle tissue
WB Recommended dilution:	1000-5000



