

## GIMAP2 抗原（重组蛋白）

中文名称：GIMAP2 抗原（重组蛋白）

英文名称：GIMAP2 Antigen (Recombinant Protein)

别名：IAN12; IMAP2; HIMAP2

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-200 amino acids of human GIMAP2

技术规格：

<b>Full name:</b>	GTPase, IMAP family member 2
<b>Synonyms:</b>	IAN12; IMAP2; HIMAP2
<b>Swissprot:</b>	Q9UG22
<b>Gene Accession:</b>	BC013934
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a protein belonging to the GTP-binding superfamily and to the immuno-associated nucleotide (IAN) subfamily of nucleotide-binding proteins. In humans, the IAN subfamily genes are located in a cluster at 7q36.1.