

## 兔抗 FIBIN 多克隆抗体

中文名称：兔抗 FIBIN 多克隆抗体

英文名称：Anti-FIBIN rabbit polyclonal antibody

别名：fin bud initiation factor homolog

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FIBIN

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

|                                    |   |
|------------------------------------|---|
| <b>Background:</b>                 | The function of this protein remains unknown. |
| <b>Applications:</b>               | ELISA, IHC                                    |
| <b>Name of antibody:</b>           | FIBIN   |
| <b>Immunogen:</b>                  | Fusion protein of human FIBIN                 |
| <b>Full name:</b>                  | fin bud initiation factor homolog             |
| <b>SwissProt:</b>                  | Q8TAL6  |
| <b>ELISA Recommended dilution:</b> | 5000-10000                                    |
| <b>IHC positive control:</b>       | Human ovarian cancer                          |
| <b>IHC Recommend dilution:</b>     | 150-300                                       |

