

## 兔抗 FIS1 多克隆抗体

中文名称：兔抗 FIS1 多克隆抗体

英文名称：Anti-FIS1 rabbit polyclonal antibody

别名：TTC11; CGI-135

相关类别：一抗

储存：冷冻 (-20℃)

抗原：FIS1

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The balance between fission and fusion regulates the morphology of mitochondria. TTC11 is a component of a mitochondrial complex that promotes mitochondrial fission (James et al., 2003 [PubMed 12783892]).
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	FIS1
<b>Immunogen:</b>	Fusion protein of human FIS1
<b>Full name:</b>	fission, mitochondrial 1
<b>Synonyms:</b>	TTC11; CGI-135
<b>SwissProt:</b>	Q9Y3D6
<b>ELISA Recommended dilution:</b>	2000-5000

<b>IHC positive control:</b>	Human liver cancer and Human brain
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	17 kDa
<b>WB Positive control:</b>	Raji cells, mouse liver and heart tissue
<b>WB Recommended dilution:</b>	500-2000



