

兔抗 FIZ1 多克隆抗体

中文名称：兔抗 FIZ1 多克隆抗体

英文名称：Anti-FIZ1 rabbit polyclonal antibody

别名：FLT3 interacting zinc finger 1; ZNF798

抗原：FIZ1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes zinc finger protein, which interacts with a receptor tyrosine kinase involved in the regulation of hematopoietic and lymphoid cells. This gene product also interacts with a transcription factor that regulates the expression of rod-specific genes in retina.
Applications:	ELISA, IHC
Name of antibody:	FIZ1
Immunogen:	Synthetic peptide of human FIZ1
Full name:	FLT3 interacting zinc finger 1
Synonyms:	ZNF798
SwissProt:	Q96SL8

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human liver cancer
IHC Recommend dilution:	40-200



