

GOLGA7 抗原（重组蛋白）

中文名称：GOLGA7 抗原（重组蛋白）

英文名称：GOLGA7 Antigen (Recombinant Protein)

别名：GCP16; GOLGA7A; HSPC041; GOLGA3AP1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-137 amino acids of human GOLGA7

技术规格：

Full name:	golgin A7
Synonyms:	GCP16; GOLGA7A; HSPC041; GOLGA3AP1
Swissprot:	Q7Z5G4
Gene Accession:	BC012032
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May be involved in protein transport from Golgi to cell surface. The ZD HHC9-GOLGA7 complex is a palmitoyltransferase specific for HRAS and NRAS.