

GPKOW 抗原（重组蛋白）

中文名称：GPKOW 抗原（重组蛋白）

英文名称：GPKOW Antigen (Recombinant Protein)

别名：T54; Spp2; GPATC5; GPATCH5

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 277-476 amino acids of human GPKOW

技术规格：

Full name:	G-patch domain and KOW motifs
Synonyms:	T54; Spp2; GPATC5; GPATCH5
Swissprot:	Q92917
Gene Accession:	BC000397
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a putative RNA-binding protein containing G-patch and KOW (Kyprides, Ouzounis, Woese) domains. The encoded protein interacts directly with protein kinase A and protein kinase X and is also found associated with the spliceosome.