

GSG1L 抗原（重组蛋白）

中文名称：GSG1L 抗原（重组蛋白）

英文名称：GSG1L Antigen (Recombinant Protein)

别名：PRO19651

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 156-331 amino acids of human GSG1L

技术规格：

Full name:	GSG1 like
Synonyms:	PRO19651
Swissprot:	Q6UXU4
Gene Accession:	BC016460
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	As a component of the inner core of AMPAR complex, modifies AMPA receptor (AMPA) gating.